

Result

1 Project information

The basic information of the project is shown below:

- Project ID: F20FTSEUHT0946-01_HOMaltyR
- Product name: filter
- Sample size: 8
- Library type: DNBSEQ Eukaryotic Strand-specific mRNA library
- Sequencing Platform: DNBseq
- Sequencing read Length: PE100
- Clean fastq phred quality score encoding:Phred+33

2 Data production

After sequencing, the raw reads were filtered. Data filtering includes removing adaptor sequences, contamination and low-quality reads from raw reads .

Table 1 Statistics of clean data					
Sample Name	Clean Reads	Clean Base	Read Length	Q20(%)	GC(%)
1	24,063,106	4,812,621,200	PE100	97.58	50.74
3	24,055,454	4,811,090,800	PE100	97.50	50.99
45_ShAP	24,056,328	4,811,265,600	PE100	98.28	53.65
5	24,070,060	4,814,012,000	PE100	97.59	50.94
7	24,048,479	4,809,695,800	PE100	97.67	51.60
K20	24,053,822	4,810,764,400	PE100	97.55	53.23
N1	24,032,872	4,806,574,400	PE100	98.33	51.68
R6	24,024,971	4,804,994,200	PE100	98.36	52.07

- Sample Name : Sample Name
- Clean Reads : Clean reads
- Clean Base : Clean bases
- Read Length : Read length
- Q20(%) : Proportion of Q20
- GC(%) : Proportion of GC

3 Quality control

The quality of data was examined after filtering.

3.1 The distribution of base percentage and qualities along reads after data filtering

































4 References

[1] SOAPnuke: a MapReduce acceleration-supported software for integrated quality control and preprocessing of high-throughput sequencing data.

Chen Y, Chen Y, Shi C, et al.

PMID: 29220494 PMCID: PMC5788068 DOI: 10.1093/gigascience/gix120



Methods

1 Experimental procedure

The library construction method and sequencing process are carried out according to the following steps:



- 1. Take a certain amount of total RNA samples, and use oligo dT beads to enrich mRNA with poly A tail;
- 2. mRNA molecules were fragmented into small pieces;
- 3. The fragmented mRNA was synthesized into first strand cDNA using random primers;
- 4. The second strand cDNA was synthesized with dUTP instead of dTTP;
- 5. The synthesized cDNA was subjected to end-repair and 3' adenylated. Adaptors were ligated to the ends of these 3' adenylated cDNA fragments;
- 6. Digest the U-labeled second-strand template with Uracil-DNA-Glycosylase (UDG) and perform PCR amplification;
- 7. Library quality control;
- 8. Library circularization;
- 9. The library was amplified to make DNA nanoball (DNB);
- 10. Sequencing on DNBSEQ (DNBSEQ Technology) platform.

2 Bioinformatic analysis workflow



Figure 2 Bioinformatic analysis workflow

3 Parameters for data filtering

raw data with adapter sequences or low-quality sequences was filtered. We first went through a series of data processing to remove contamination and obtain valid data. This step was completed by SOAPnuke software developed by BGI.

SOAPnuke software filter parameters: " -n 0.01 -l 20 -q 0.4 --adaMR 0.25 --ada_trim --polyX 50", steps of filtering: SOAPnuke software filtering parameters: '-n 0.01 -l 20\\-q 0.4 --adaMR 0.25 -- ada_trim --polyX 50', filtering steps:

- 1. Filter adapter: $\ \$ if the sequencing read matches 25.0% or more of the adapter sequence (maximum $\ \ \ 2$ base mismatches are allowed), cut the adapter;
- Filter low-quality data:\\ if the bases with a quality value of less than 20 in the sequencing read account\ \ for 40.0% or more of the entire read, delete the entire read;
- 3. Remove N: $\ \$ if the N content in the sequencing read accounts for 1.0% or more of the entire $\ \$ read, delete the entire read;
- 4. Remove polyX: if the length of polyX (X $\ C$ or C) in the sequencing read exceeds 50bp, delete the entire $\ \$

4 References

[1] SOAPnuke: a MapReduce acceleration-supported software for integrated quality control and preprocessing of high-throughput sequencing data.

Chen Y, Chen Y, Shi C, et al.

PMID: 29220494 PMCID: PMC5788068 DOI: 10.1093/gigascience/gix120



Help

1 FASTQ format description

Images generated by sequencers are converted by base calling into nucleotide sequences, which are called raw data or raw reads and are stored in FASTQ format. FASTQ files are text files that store both reads sequences and their corresponding quality scores. Each read is described in four lines as follows:

The first line is the sequence identifier and related description information, starting with'@'; the second line is the base sequence information; the third line starts with'+', followed by the sequence identifier, description information, or nothing; The four lines are quality information, which corresponds to the sequence in the second line. Each base has a quality score. Depending on the scoring system, each character represents a different quality value.

The figure below shows the concise correspondence between the sequencing error rate and the sequencing quality value. Specifically, if the sequencing error rate is denoted by E and the base quality value is denoted by SQ, there are the following relationships:

$$SQ = -10 \times (\log \frac{E}{1-E}) / (\log 10)$$
$$E = \frac{Y}{1+Y}$$
$$Y = \frac{SQ}{e^{-10 \times \log 10}}$$

- For the quality system data with a sequencing quality value of 33: the sequencing quality value of the base = the ASCII value corresponding to the quality information character -33, for example, the ASCII value corresponding to A is 65, then the corresponding base quality value is 65-33 = 32. The base quality value of the DNBSEQ sequencing platform ranges from 2 to 42.
- 2. For quality system data with a sequencing quality value of 64: the sequencing quality value of the base = the ASCII value corresponding to the quality information character -64, for example, the corresponding ASCII value of c is 99, then the corresponding base quality value is 99-64 =35. The base quality value of the DNBSEQ sequencing platform ranges from 2 to 43.

Table 1 A summary table of the relationship between sequencing error rate and sequencing quality						
Sequencing error rate	Sequencing quality value	Character(Phred64)	Character(Phred33)			
5%	13	Μ				
1%	20	Т	5			
0.1%	30	٨	?			